


# NurseCare interface for connection of wired and wireless units with RFID function

**Art.: ERTW0015 (BTID)**



## Description:

NurseCare interface for connection of wired and wireless units. It's used like a bridge between Nursetab terminal and other wired/wireless NurseCare devices. Includes interface for connection of wired and wireless units with RFID antenna.

## Technical specification:

- Includes connection to room terminal
- Receives signal from NurseCare units for staff protection
- Receives signal from NurseCare nurse call system units
- Integrated RFID reader with RFID antenna for identification of staff
- Includes location functionality
- Numbers of Inputs: 6

- 
- Numbers of outputs: 6
  - Dimensions (W: 80mm / H: 40mm / D: 12mm)
  - Weight: 12g
  - Frequency: 2.4 GHz ISM (2.402 GHz -2.480 GHz Utilized)
  - Transmission power 10 mW
  - Voltage IN: 24VDC
  - Voltage OUT: 24VDC
  - Operating temperature: +5°C - to +55°C
  - Operating humidity: Up to 85%, non-condensing
  - Supported RFID card types: NFC Forum Tag Type 1, 2 and 4 (ISO/IEC 14443A)
  - Compliant with CE certification requirements
  - Compliant with DIN VDE 0834 certification requirements